



HOMOGENEOUS ENZYMIC HOMOCYSTEINE REAGENT

Hitachi 917 Procedure Parameters:

Test: HCY* Type: Ser/PI

ANALYZE

Assay time/Point	[2 Point End]/[10]/[19] [34] [0] [0]
Wave (2nd/Primary)	[376]/[340]
S. Vol (Normal)	[18.6]
Reagent (R1) T1	[270] [0] [000000]
Reagent (R2) T2	[33] [0] [000000]
Reagent (R3) T3	[27] [0] [000000]
Abs. Limit	[32000] [Decrease] 2Tests
Prozone Limit	[-32000] [0] [Lower]
Cell Detergent	[Detergent 1]

CALIB

Calibration Type	[Linear]
Point	[2 Span Point]
Weight	[0]
Auto Calibration	
Full	[120]
SD Limit	[100]
Duplicate Limit	[10%] [32000 Abs]
Sensitivity Limit	[-99999] [99999]
S1 Abs limit	[-32000] [32000]

RANGE

Application Code*	[]	Unit [μmol/L]
Control Interval*	[]	
Instrument Factor	(Y=aX+b) a=[1.0] b= [0.0]	
Technical Limit	[1.0] [50.0]	
Repeat Limit*	[-99999] [99999]	

OTHERS

<Standard>	(1)	(2)
Calibration Code*	[]	[]
Concentration**	[0.00]	[**]
Position*	[]	[]
Sample Volume	[18.6]	[18.6]

RESULTS

Results are printed out by the Hitachi 917 in μmol/L.

* User Defined

** Enter Values on Calibrator Vials

Note: Samples with values greater than 50 μmol/L should be diluted 1:3 with saline and rerun. Multiply results x3.