



HOMOGENEOUS ENZYMIC HOMOCYSTEINE REAGENT

Roche Modular System Procedure Parameters:

Test No [*]	Test Name [HCY]	Test Code [*]
Sample Type	Ser/PI	
Assay/ time/Point	[2 Point End]/[10]/[19] [34] [0] [0]	
Wave (Primary/2nd)	[340]/[376]	
S. Vol (Normal/Dec/Inc)	[18.6] / [] / []	
Reagent (R1)	[270] [0] [000000]	
Reagent (R2)	[33] [0] [000000]	
Reagent (R3)	[27] [0] [000000]	
Abs. Limit	[32000] [Decrease]	
Prozone Limit	[-32000] [Lower]	
Point	[]	
Tech Limit (Low/High)	[1.0]/[50.0]	
Rpt Limit (Low/High)	[-99999]/[999999]	
Calibration Method	[Linear]	
Point	[2]	
Span Point	[2]	
Update Type		
Duplicate Limit (%)	[10]	
Dup Limit (ABS)	[32000]	
SD Limit	[100]	
Sensitivity Limit	[-99999] [999999]	
S1 Abs limit	[-32000] [32000]	
<Standard>	(1)	(2)
Concentration	[0.0]	[**]
Automatic Rerun	Yes	
Auto Masking	Yes	

* User Defined

** Enter Values on Calibrator Vials

Note: Samples with values greater than 50 µmol/L should be diluted 1:3 with the saline and rerun.
Multiply results by 3.

RESULTS:

Results are printed out by the Roche Modular in µmol/L.